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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

### **APPLICATION AND DESCRIPTION OF THE DEVICE**

Professional microscope with magnification up to 1000x and LED backlight allows you to see the smallest scratch or fingerprint on your hand. It allows you to see much more than is possible with the naked eye. It uses a USB interface to connect to a computer or laptop, thanks to which it provides a continuous preview on the screen and has the ability to record video or take photos. A digital microscope is an invaluable tool that will easily make any training, biology lesson or presentation more attractive. Thanks to the ability to transfer a digital image in real time to a computer, and from there, for example, to a wall projector, our presentations will enjoy great interest. Compact dimensions allow you to take it anywhere. It will find use in collecting, electronics, jewelry or telephone services. The microscope is ideal for all precise work requiring visual focus. The high quality of the lenses used in the production of sensors makes the image sharp and contrasting. Thanks to 8 strong LEDs and adjustable light intensity, observations are possible both day and night!

For safety and CE certification reasons, the product may not be rebuilt or modified in any way. If the product is used for purposes other than those previously described, the product may be damaged. Improper use may also cause hazards such as short circuits, fire, electric shock, etc.

### **COMMISSIONING/INSTALLATION INSTRUCTIONS**

- Note: Packaging material must be completely removed from the product.
- Make sure the components are assembled correctly.
- Cover sensitive surfaces before placing the product on the table.
- Place the product near a wall for best results. Note: The product is not intended for ceiling projection.
- Connect the USB cable to a USB power source and connect the DC plug of the USB cable to the power input port of the product.

### **TIP**

- The device must be checked for completeness of delivery and any visible damage.
- In the event of an incomplete delivery or if you notice damage due to faulty packaging or transport, please contact the service hotline.

### **TECHNICAL DATA**

- Image sensor: HD color CMOS
- Controller: High Speed DSP
- Maximum sharpness: up to 8mm
- Photos: via software or by button on the device (SNAP)
- Video format: AVI 30fps
- Photo and video resolution: standard (640x480), maximum (1600x1200)
- Matrix: 2Mpix
- Snap Shot software and hardware included
- Lighting: 8 x SMD white LEDs with adjustable light intensity
- Digital zoom: 6x - sequential mode
- Optical zoom: up to 1000x adjustable by knob
- Power supply: USB port 5V
- Microscope dimensions: 110mm x 34mm
- USB cable length: 135cm

### **SAFETY INSTRUCTIONS**

- This product may be used by children aged 8 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- Never immerse the device in water.
- Use a soft, damp cloth for cleaning. Avoid using harsh cleaning agents that may damage the device.
- Do not use a damaged device.
- The product is not intended for use by children.
- Make sure that packaging materials are not left unattended. Children may start playing with them, which is dangerous.
- Protect the product from extreme temperatures, direct sunlight, strong vibrations, high humidity, moisture, flammable gases, vapors and solvents.
- Do not expose the product to mechanical stress.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - does not function properly, - has been stored for an extended period under unfavourable conditions or - has been excessively stressed during transport.

- It is forbidden to use the product if any of its parts are damaged. In the event of damage to the cable, it is forbidden to make repairs yourself.
- Make sure the LED light is turned off when not in use to save batteries.
- Do not point the LED light directly into the eyes to avoid eye irritation.
- Avoid touching the lens surface with your fingers to avoid leaving streaks. Clean them with a soft microfiber cloth.
- This product contains small parts that may present a choking hazard. Keep out of reach of children.
- When the magnifier is not used for a long time, remove the batteries to prevent them from leaking.
- Do not disassemble the device yourself.
- Do not use old and new batteries, different brands or chemicals.
- Do not charge, damage or open the battery.
- Keep batteries out of the reach of children to prevent accidental swallowing.
- Do not swallow the battery. If the battery is swallowed, consult a doctor immediately.
- Always use the product as intended.

### **BATTERIES/ACCUMULATORS SAFETY INSTRUCTIONS**

- Batteries/rechargeable batteries should be kept out of reach of children. If swallowed, consult a doctor immediately!
- Disposable batteries must not be recharged. Batteries/accumulators must not be short-circuited and/or opened. This may cause overheating, fire or explosion.
- Never throw batteries/accumulators into fire or water.
- Batteries/accumulators should never be subjected to mechanical stress.
- Risk of acid leakage from batteries/accumulators.
- Avoid extreme conditions and temperatures that may affect batteries / accumulators, e.g. radiators / direct sunlight.
- If batteries/accumulators leak, avoid contact of skin, eyes and mucous membranes with chemicals! Immediately rinse affected areas with clean water and seek medical attention!
- Leaking or damaged batteries/rechargeable batteries can cause chemical burns when in contact with skin. Therefore, suitable protective gloves should be worn in such cases.

#### **TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING**

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

#### **DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT**

Used electrical and electronic products, due to environmental protection, must not be thrown into household waste with regular municipal waste, but disposed of properly. Information on collection points and their opening hours is provided by the relevant office.

This product complies with the requirements of the relevant European and national directives that apply to it. The product meets European and national requirements on the safety of devices and products.

This product complies with the requirements of the relevant European and national directives that apply to it, RoHS.

Do not look into a light source when it is on. It may be harmful to your eyes.

#### **Improper disposal of batteries/accumulators poses a risk to the environment!**

Batteries/accumulators should not be disposed of with household waste. They may contain harmful heavy metals and should be treated as special waste. The chemical symbols for heavy metals are as follows: Cd = cadmium, Hg = mercury, Pb = lead.

Therefore, used batteries / accumulators should be taken to municipal hazardous waste collection points.

The built-in battery must not be disassembled for disposal. The product must be handed over in its entirety to a collection point for used electronics.

We reserve the right to make changes to the text, design and technical data of the product without notice.

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