- (PL) INSTRUKCJA OBSŁUGI - WSKAZÓWKI DOTYCZĄCE OBSŁUGI I BEZPIECZEŃSTWA**
TERMOMETR, SYMBOL: 01631 EAN/GTIN: 5907451329719
- (EN) USER MANUAL - OPERATING AND SAFETY INSTRUCTIONS**
THERMOMETER, SYMBOL: 01631 EAN/GTIN: 5907451329719
- (DE) BEDIENUNGSANLEITUNG - BEDIENUNGS- UND SICHERHEITSHINWEISE**
THERMOMETER, SYMBOL: 01631 EAN/GTIN: 5907451329719
- (CZ) NÁVOD K OBSLUZE - PROVOZNÍ A BEZPEČNOSTNÍ RADY**
TEPLOMĚR, SYMBOL: 01631 EAN/GTIN: 5907451329719
- (FR) MANUEL D'UTILISATION - CONSEILS D'UTILISATION ET DE SÉCURITÉ**
THERMOMÈTRE, SYMBOLE : 01631 EAN/GTIN: 5907451329719
- (IT) MANUALE OPERATIVO - CONSIGLI OPERATIVI E DI SICUREZZA**
TERMOMETRO, SIMBOLO: 01631 EAN/GTIN: 5907451329719
- (ES) MANUAL DE OPERACIÓN - CONSEJOS DE OPERACIÓN Y SEGURIDAD**
TERMÓMETRO, SÍMBOLO: 01631 EAN/GTIN: 5907451329719
- (NL) BEDIENINGSHANDLEIDING - BEDIENINGS- EN VEILIGHEIDTIPS**
THERMOMETER, SYMBOL: 01631 EAN/GTIN: 5907451329719
- (SE) BRUKSANVISNING - DRIFT- OCH SÄKERHETSTIPS**
TERMOMETR, SYMBOL: 01631 EAN/GTIN: 5907451329719
- (GR) ΕΓΧΕΙΡΙΔΙΟ ΛΕΙΤΟΥΡΓΙΑΣ - ΣΥΜΒΟΥΛΕΣ ΛΕΙΤΟΥΡΓΙΑΣ ΚΑΙ ΑΣΦΑΛΕΙΑΣ**
ΘΕΡΜΟΜΕΤΡΟ, ΣΥΜΒΟΛΟ: 01631 EAN/GTIN: 5907451329719
- (RO) MANUAL DE OPERARE - SFATURI DE OPERAREA SI SIGURANTA**
TERMOMETRU, SIMBOL: 01631 EAN/GTIN: 5907451329719
- (PT) MANUAL DE OPERAÇÃO - DICAS DE OPERAÇÃO E SEGURANÇA**
TERMÔMETRO, SÍMBOLO: 01631 EAN/GTIN: 5907451329719
- (BG) РЪКОВОДСТВО ЗА ЕКСПЛОАТАЦИЯ - СЪВЕТИ ЗА ЕКСПЛОАТАЦИЯ И БЕЗОПАСНОСТ**
ТЕРМОМЕТЪР, СÍМБОЛ: 01631 EAN/GTIN: 5907451329719
- (HU) HASZNÁLATI ÚTMUTATÓ – HASZNÁLATI ÉS BIZTONSÁGI TIPP EK**
HŐMÉRŐ, SZIMBÓLUM: 01631 EAN/GTIN: 5907451329719
- (DK) BETJENINGSVEJLEDNING - BETJENINGS- OG SIKKERHEDSTIPS**
TERMOMETR, SYMBOL: 01631 EAN/GTIN: 5907451329719
- (SK) NÁVOD NA OBSLUHU - RADY NA PREVÁDZKU A BEZPEČNOSŤ**
TEPLOMER, SYMBOL: 01631 EAN/GTIN: 5907451329719
- (FI) KÄYTTÖOHJE - KÄYTTÖ- JA TURVALLISUUSVINKKEJÄ**
LÄMPÖMITTARI, SYMBOLO: 01631 EAN/GTIN: 5907451329719
- (LT) NAUDOJIMO VADOVAS – NAUDOJIMO IR SAUGOS PATARIMAI**
TERMOMETRAS, SIMBOLIS: 01631 EAN/GTIN: 5907451329719
- (LV) LIETOŠANAS ROKASGRĀMATA – LIETOŠANAS UN DROŠĪBAS PADOMI**
TERMOMETRS, SIMBOLS: 01631 EAN/GTIN: 5907451329719
- (EE) KASUTUSJUHEND – KASUTUS- JA OHUTUSJUHEND**
TERMOMEETER, SÜMBOOL: 01631 EAN/GTIN: 5907451329719
- (SI) PRIROČNIK ZA UPORABO - NASVETI ZA UPORABO IN VARNOST**
TERMOMETR, SIMBOL: 01631 EAN/GTIN: 5907451329719
- (IE) LÁMHLEABHAR OIBRIÚCHÁIN - TIPS OIBRIÚCHÁIN AGUS SÁBHÁILTEACHTA**
TEIRMIMÉADAR, SIMBOL: 01631 EAN/GTIN: 5907451329719
- (MT) MANWAL OPERAT - OPERAZZJONI U GHAJR TAS-SIGURTA**
TERMOMETRU, SIMBOLU: 01631 EAN/GTIN: 5907451329719
- (HR) PRIRUČNIK ZA UPORABU - SAVJETI ZA UPOTREBU I SIGURNOST**
TERMOMETAR, SIMBOL: 01631 EAN/GTIN: 5907451329719
- (RU) РУКОВОДСТВО ПО ЭКСПЛУАТАЦИИ – СОВЕТЫ ПО ЭКСПЛУАТАЦИИ И БЕЗОПАСНОСТИ**
ТЕРМОМЕТР, СИМБОЛ: 01631 EAN/GTIN: 5907451329719

Importer / Importeur / Importateur / Importatore / Importador / Importör / Εισαγωγέας / Importator / Dovožce / Importador / Вносител / Importör / Importør / Dovožca / Maahantuoja / Importuotojas / Importētājs / Maaletooja / Uvoznik / Allmhaireoir / Importatur / Uvoznik / Импортёр:

VERK GROUP SIKORSKI SP.K.

WYGODY 16

05-090 PODOLSZYN NOWY

POLAND

info@verkgroup.pl

+48 888055550 / +48 22 670 90 10

internetowa-hurtownia.pl

verk.store

verk.pl

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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

APPLICATION AND DESCRIPTION OF THE DEVICE

A non-contact thermometer is an essential device in every home. It offers measurement of body temperature, ambient temperature and various objects. The thermometer measures temperature in one second. This is very convenient when measuring small, impatient children. Simply point the thermometer at the forehead from a distance of three to five centimetres to get the result instantly. The thermometer is equipped with a large, legible LCD display and white backlight, which makes reading it much easier when it is dark. The correct temperature is signalled with one sound signal and fever with several short signals. The device is equipped with a memory of the last 10 measurements.

For safety and CE certification reasons, the product may not be rebuilt or modified in any way. If the product is used for purposes other than those previously described, the product may be damaged. Improper use may also cause hazards such as short circuits, fire, electric shock, etc.

Battery replacement instructions:

- Please insert the batteries correctly, paying attention to the +/- polarity.
- Only use one type of battery.
- Always replace the entire set of batteries.
- Exhausted batteries must be removed from the device and never thrown into fire.
- Never attempt to charge regular batteries.
- Store unused batteries in their packaging away from metal objects.
- Do not leave children unsupervised.

TIP

- The device must be checked for completeness of delivery and any visible damage.
- In the event of an incomplete delivery or if you notice damage due to faulty packaging or transport, please contact the service hotline.

TECHNICAL DATA

- Touchless
- Quick measurement
- Sound signal
- Hygienic measurement
- Measurement range: 3-5cm
- Measurement in Celsius or Fahrenheit
- Temperature measurement range: forehead measurement mode: 32 ~ 42.5 °C
- Forehead measurement accuracy: +/- 0.2 °C
- Ambient temperature measurement range: 0 °C ~ 100 °C
- Ambient temperature measurement accuracy: +/- 1 °C
- Number of body measurements stored: 32
- Number of saved ambient measurements: 10
- Large backlit LCD display
- Handy
- Automatic shutdown after about 15 seconds
- Power supply: 2 x AA (not included)
- Dimensions: 14.9cm x 9.6cm x 4.3cm

SAFETY INSTRUCTIONS

- This product may be used by children aged 8 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- Never immerse the device in water.
- The infrared thermometer is not waterproof. Do not use it near water or in places exposed to moisture to avoid damage to the device.
- Do not point the laser beam at people, animals, or reflective surfaces such as mirrors or metal. This may cause eye damage.
- Avoid measuring the temperature on very hot or flammable surfaces that may cause overheating or fire.
- Turn your device off after use to avoid unnecessary battery consumption.
- Keep the device out of reach of children. The laser and thermometer components may be dangerous if misused.
- Check the batteries regularly and replace them when the device stops working properly. Dispose of used batteries according to local regulations.
- Use a soft, damp cloth for cleaning. Avoid using harsh cleaning agents that may damage the device.

- Do not use a damaged device.
- The product is not intended for use by children.
- Make sure that packaging materials are not left unattended. Children may start playing with them, which is dangerous.
- Protect the product from extreme temperatures, direct sunlight, strong vibrations, high humidity, moisture, flammable gases, vapors and solvents.
- Do not expose the product to mechanical stress.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - does not function properly, - has been stored for an extended period under unfavourable conditions or - has been excessively stressed during transport.
- It is forbidden to use the product if any of its parts are damaged. In the event of damage to the cable, it is forbidden to make repairs yourself.
- This product contains small parts that may present a choking hazard. Keep out of reach of children.
- Do not disassemble the device yourself.
- Do not use old and new batteries, different brands or chemicals.
- Do not charge, damage or open the battery.
- Keep batteries out of the reach of children to prevent accidental swallowing.
- Do not swallow the battery. If the battery is swallowed, consult a doctor immediately.
- Always use the product as intended.

BATTERIES/ACCUMULATORS SAFETY INSTRUCTIONS

- Batteries/rechargeable batteries should be kept out of reach of children. If swallowed, consult a doctor immediately!
- Disposable batteries must not be recharged. Batteries/accumulators must not be short-circuited and/or opened. This may cause overheating, fire or explosion.
- Never throw batteries/accumulators into fire or water.
- Batteries/accumulators should never be subjected to mechanical stress.
- Risk of acid leakage from batteries/accumulators.
- Avoid extreme conditions and temperatures that may affect batteries / accumulators, e.g. radiators / direct sunlight.
- If batteries/accumulators leak, avoid contact of skin, eyes and mucous membranes with chemicals! Immediately rinse affected areas with clean water and seek medical attention!
- Leaking or damaged batteries/rechargeable batteries can cause chemical burns when in contact with skin. Therefore, suitable protective gloves should be worn in such cases.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT

Used electrical and electronic products, due to environmental protection, must not be thrown into household waste with regular municipal waste, but disposed of properly. Information on collection points and their opening hours is provided by the appropriate office.

This product complies with the requirements of the relevant European and national directives that apply to it. The product meets European and national requirements on the safety of devices and products.

This product complies with the requirements of the relevant European and national directives that apply to it, RoHS.

Improper disposal of batteries/accumulators poses a risk to the environment!

Batteries/accumulators should not be disposed of with household waste. They may contain harmful heavy metals and should be treated as special waste. The chemical symbols for heavy metals are as follows: Cd = cadmium, Hg = mercury, Pb = lead.

Therefore, used batteries / accumulators should be taken to municipal hazardous waste collection points.

The built-in battery must not be disassembled for disposal. The product must be handed over in its entirety to a collection point for used electronics.

We reserve the right to make changes to the text, design and technical data of the product without notice.

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