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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

APPLICATION AND DESCRIPTION

The acupressure mat is a rectangular pad equipped with hundreds of spikes arranged in a specific way. Each spike creates micro-pressure points that stimulate receptors on the skin. The mat has 230 massaging rosettes and a large surface area, thanks to which it has increased therapeutic effectiveness. The mattress is covered with a cotton cover fastened with a comfortable Velcro, which allows you to keep the mat clean. The spikes placed in the mattress, under the influence of the body weight, exert pressure on the skin. This increases the blood flow to the pressed area of the body, relaxes the muscles and improves vital functions.

TIP

- The device must be checked for completeness of delivery and any visible damage.
- In the event of an incomplete delivery or if you notice damage due to faulty packaging or transport, please contact the service hotline.

TECHNICAL DATA

- 230 massage rosettes
- Diameter of the rosette: 3cm
- Number of spikes on the rosette: 27
- Removable cover
- Length: 65cm
- Width: 41cm
- Height: 2cm
-

SAFETY INSTRUCTIONS

- Start with short sessions (5-10 minutes) to allow your body to get used to the spikes.
- You can use a thin T-shirt or towel to reduce the intensity of the pressure.

Avoid use in case of:

- Open wounds, skin lesions or infections.
- Skin diseases or hypersensitivity to touch.
- Circulatory problems, haemophilia or use of anticoagulant medications.
- Pregnancy (consultation with a doctor is recommended).

During use:

- Make sure the mat is placed on a flat surface.
- Avoid sudden pressure – relax gradually.
- Do not exceed 20-30 minutes in one session.

After use:

- Get up slowly to avoid dizziness.
- Clean the mat by wiping it with a damp cloth or according to the manufacturer's instructions.

Comments :

- Regular use of the mat can bring the best results, but always listen to your body's signals.
- If you have any doubts or health problems, consult your doctor before use.
- Always use the product as intended.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

The product meets European and national requirements on the safety of devices and products.

We reserve the right to make changes to text, design and product data without notice.

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