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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

APPLICATION AND DESCRIPTION OF THE DEVICE

This product is intended for home use. This product is not intended for commercial use or other applications.

The muscle electrostimulator is a therapy device designed for home use. The electrostimulation device can be used in sports, rehabilitation, body care and for beauty improvement.

It is effective in reducing fat tissue or eliminating cellulite. It also helps in treating pain. Hence, muscle electrostimulation is an extremely useful procedure in the context of convalescence after many injuries. Athletes also appreciate the devices due to better regeneration and elimination of pain and soreness, which were very bothersome the day after the training.

The set includes a pair of flip-flops that act as foot electrodes, so you can massage your tired feet after an exhausting day at work.

Muscle training and increasing their endurance

Relaxation

Muscle fatigue

Muscle regeneration after exercise

Back pain

Joint pain

Neuralgia

For post-traumatic pain of the musculoskeletal system

In case of blood circulation disorders

In chronic pain conditions of various causes

For safety and CE certification reasons, the product may not be rebuilt or modified in any way. If the product is used for purposes other than those previously described, the product may be damaged. Improper use may also cause hazards such as short circuits, fire, electric shock, etc.

STARTING/ USING FUNCTION BUTTONS

- Note: Packaging material must be completely removed from the product.
- Make sure the components are assembled correctly.
- Electrodes should be stuck to the body in selected places:

SHOULDER button - massages the shoulders and back

WAIST button - massages the waist and surrounding areas

JOINT button - stimulates the joints

HAND&FOOT button - massages feet and hands

LEG button - legs

Using the buttons: **PRESS, RUB, VIBRATE, THUMB**, we change the type of massage performed.

Inserting / Replacing the battery

- Open the battery compartment cover and remove it.
- Remove the used batteries. Insert new batteries into the battery compartment.
- Pay attention to the orientation of the poles (+ and -).
- Close the battery compartment cover.

Connecting to USB Power

- Remove the batteries from the battery compartment before establishing a USB connection.
- Connect the USB cable to a USB power source.
- Connect the DC plug of the USB cable to the power socket of the product.

Selecting the electrical power supply

If the DC plug is not connected to the product, power is supplied by batteries.

Battery replacement instructions:

- Please insert the batteries correctly, paying attention to the +/- polarity.
- Only use one type of battery.
- Always replace the entire set of batteries.
- Exhausted batteries must be removed from the device and never thrown into fire.
- Never attempt to charge regular batteries.
- Store unused batteries in their packaging away from metal objects.
- **Do not leave children unsupervised.**

TIP

- The device must be checked for completeness of delivery and any visible damage.
- In the event of an incomplete delivery or if you notice damage due to faulty packaging or transport, please contact the service hotline.

TECHNICAL DATA

- Product code: 15833
- Easy-to-read LCD display
- Possibility to connect 4 electrodes or a pair of flaps
- Velcro flaps: insert length 26cm
- Speed and intensity adjustment
- Intuitive operation
- Effective selection of programs
- Selecting the massage type1x connecting cable (4 CLIPS)
- 1x connecting wire (2 CLIPS)
- Electrode pad
- USB charger and power cable
- Battery powered: 3x AA (not included)
- The product is CE certified
- Note! The device should not be used by people with a pacemaker, insulin pump or pregnant women.

SAFETY OF USE

Before determining the place of application of electrodes, always remember the general rules that do not allow electrodes to be placed in the following places:

On opposite sides of the braincase in such a way that the electrical impulse runs through the brain

On opposite sides of the chest so that the impulse runs through the mediastinum

Directly on the wounds

Directly above the carotid sinus (anterior lateral middle part of the neck)

On the eyelids and above the eye socket

Directly to the sites of inflammation or infection

Inside the oral cavity on the mucosa

On densely hairy skin (skin can be shaved)

Note! The device should not be used by people with a pacemaker, insulin pump or pregnant women.

SAFETY INSTRUCTIONS

- This product may be used by children aged 8 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- Children are not allowed to play with the product.
- Cleaning and maintenance may not be performed by children without supervision.
- Never immerse the device in water.
- For cleaning, use a damp cloth or a mild detergent.
- Do not use a damaged device.
- Make sure that packaging materials are not left unattended. Children may start playing with them, which is dangerous.
- Protect the product from extreme temperatures, direct sunlight, strong vibrations, high humidity, moisture, flammable gases, vapors and solvents.
- Do not expose the product to mechanical stress.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - does not function properly, - has been stored for an extended period under unfavourable conditions or - has been excessively stressed during transport.
- It is forbidden to use the product if any of its parts are damaged. In the event of damage to the cable, it is forbidden to make repairs yourself.
- Do not disassemble the device yourself.
- Protect the product from moisture.
- Make sure not to damage the USB cable with sharp edges or hot objects. Unwind the USB cable completely before use.
- The included USB cable is only suitable for this product.
- Immediately turn off the product and remove the cable from the product if it smells burning or smoke is visible. Have the product checked by a qualified technician before using it again.
- If a USB power supply is used, the socket must always be easily accessible so that the USB power supply can be easily removed from the socket in an emergency. Please also observe the manufacturer's operating instructions.
- Always use the product as intended.

BATTERIES/ACCUMULATORS SAFETY INSTRUCTIONS

- Batteries/rechargeable batteries should be kept out of reach of children. If swallowed, consult a doctor immediately!
- Disposable batteries must not be recharged. Batteries/accumulators must not be short-circuited and/or opened. This may cause overheating, fire or explosion.
- Never throw batteries/accumulators into fire or water.
- Batteries/accumulators should never be subjected to mechanical stress.
- Risk of acid leakage from batteries/accumulators.
- Avoid extreme conditions and temperatures that may affect batteries / accumulators, e.g. radiators / direct sunlight.
- If batteries/accumulators leak, avoid contact of skin, eyes and mucous membranes with chemicals! Immediately rinse affected areas with clean water and seek medical attention!
- Leaking or damaged batteries/rechargeable batteries can cause chemical burns when in contact with skin. Therefore, suitable protective gloves should be worn in such cases.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT

Used electrical and electronic products, due to environmental protection, must not be thrown into household waste with regular municipal waste, but disposed of properly. Information on collection points and their opening hours is provided by the relevant office.

This product complies with the requirements of the relevant European and national directives that apply to it. The product meets European and national requirements on the safety of devices and products.

This product complies with the requirements of the relevant European and national directives that apply to it, RoHS.

Do not look into a light source when it is on. It may be harmful to your eyes.

Improper disposal of batteries/accumulators poses a risk to the environment!

Batteries/accumulators should not be disposed of with household waste. They may contain harmful heavy metals and should be treated as special waste. The chemical symbols for heavy metals are as follows: Cd = cadmium, Hg = mercury, Pb = lead.

Therefore, used batteries / accumulators should be taken to municipal hazardous waste collection points.

The built-in battery must not be disassembled for disposal. The product must be handed over in its entirety to a collection point for used electronics.

We reserve the right to make changes to the text, design and technical data of the product without notice.

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