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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

APPLICATION

The EMS muscle electrostimulator is a wireless training device that uses electrical impulses to activate the abdominal muscles, buttocks, thighs and arms. It is intended for both physically active people and those who want to take care of their figure without intense physical effort.

The device is intended for home use only and is not intended for medical or commercial use and is not a substitute for medical therapy. The device is intended for indoor use in dry rooms. It is not intended for medical or therapeutic use or for commercial use.

For safety and CE certification reasons, the product may not be rebuilt or modified in any way. If the product is used for purposes other than those previously described, the product may be damaged. Improper use may also cause hazards such as short circuits, fire, electric shock, etc.

ACTIVATION

- Note: Packaging material must be completely removed from the product.
- Make sure that the device components are properly assembled.
- Before first use, charge the device with a USB cable. Connect the USB cable to a USB power source and connect the USB cable plug to the USB input port of the product.
- Disconnect USB cables from devices.
- Attach the control module (black element with buttons) to the electrodes (mat with contacts).
- Remove the protective foil from the electrodes and stick the device to clean, dry skin on the selected body area (e.g. abdomen, thighs, buttocks, arms).
- Turn on the device by holding the power button (M).
- Select one of the 6 training programs using the "M" button.
- Select the operating mode and intensity using the "+" and "-" buttons. (We recommend starting with the lowest value)
- A workout typically lasts 10 minutes – the device will automatically turn off once the cycle is complete.
- Once finished, gently remove the electrodes from your skin, re-apply the protective film and store the device.

NOTE! The device should not be used by people with a pacemaker, insulin pump or pregnant women!

TIP

- The device must be checked for completeness of delivery and any visible damage.
- In the event of an incomplete delivery or if you notice damage due to faulty packaging or transport, please contact the service hotline.

TECHNICAL DATA

- 6 different programs
- 9 intensity levels
- Power supply: battery
- Charging via USB cables
- Detachable control panel
- Belt dimensions: 60cm x 20cm x 20cm
- Band dimensions: 19.5cm x 7cm
- The training lasts 10 minutes
- For exercising abdominal muscles, arms, legs

SAFETY INSTRUCTIONS

- Please read the entire instruction manual before first use.
- This product may be used by children aged 8 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- Never immerse the device in water.
- For cleaning, use a damp cloth or a mild detergent.
- Do not use a damaged device.
- The product is not intended for use by children.

- Make sure that packaging materials are not left unattended. Children may start playing with them, which is dangerous.
- Protect the product from extreme temperatures, direct sunlight, strong vibrations, high humidity, moisture, flammable gases, vapors and solvents.
- Do not expose the product to mechanical stress.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - does not function properly, - has been stored for an extended period under unfavourable conditions or - has been excessively stressed during transport.
- It is forbidden to use the product if any of its parts are damaged. In the event of damage to the cable, it is forbidden to make repairs yourself.
- Do not disassemble the device yourself.
- Protect the product from moisture.
- Make sure not to damage the USB cable with sharp edges or hot objects. Fully unwind the USB cable before use.
- The included USB cable is only suitable for this product.
- Immediately turn off the product and remove the cable from the product if it smells burning or smoke is visible. Have the product checked by a qualified technician before using it again.
- If a USB power supply is used, the socket must always be easily accessible so that the USB power supply can be easily removed from the socket in an emergency. Please also observe the manufacturer's operating instructions.
- Before each charging process, check the USB cable and device for any damage. Never use or charge the device if you find any damage.
- Before use, make sure that the existing mains voltage matches the required operating voltage of the device.
- Do not use the device if you have a pacemaker or other implanted electronic devices.
- Avoid use on wounds, irritated skin, areas with eczema or infections.
- Do not use the device while sleeping.
- Do not use while driving or operating machinery.
- Do not immerse the device in water or clean it under running water.
- Consult your doctor before first use, especially if you have heart disease, epilepsy or other neurological conditions.
- The skin should be clean, dry and free of creams, oils or sweat – these can reduce impulse conduction.
- Do not place electrodes near the heart , neck , head , face , spine , or knee or elbow joints .
- Never use the device on both sides of your chest at the same time - this may affect your heart rhythm.
- Do not stick electrodes on bones – it is best to apply them to the muscle area .
- Do not use the device while charging.
- Muscles need recovery – do not use for long periods of time without a break.
- Prolonged use may cause pain, numbness or skin irritation .
- Always use the product as intended.

BATTERIES/ACCUMULATORS SAFETY INSTRUCTIONS

- Batteries/rechargeable batteries should be kept out of reach of children. If swallowed, consult a doctor immediately!
- Disposable batteries must not be recharged. Batteries/accumulators must not be short-circuited and/or opened. This may cause overheating, fire or explosion.
- Never throw batteries/accumulators into fire or water.
- Batteries/accumulators should never be subjected to mechanical stress.
- Risk of acid leakage from batteries/accumulators.
- Avoid extreme conditions and temperatures that may affect batteries / accumulators, e.g. radiators / direct sunlight.
- If batteries/accumulators leak, avoid contact of skin, eyes and mucous membranes with chemicals! Immediately rinse affected areas with clean water and seek medical attention!
- Leaking or damaged batteries/rechargeable batteries can cause chemical burns when in contact with skin. Therefore, suitable protective gloves should be worn in such cases.
- The product has a built-in battery that cannot be replaced by the user.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT

Used electrical and electronic products, due to environmental protection, must not be thrown into household waste with regular municipal waste, but disposed of properly. Information on collection points and their opening hours is provided by the relevant office.

This product complies with the requirements of the relevant European and national directives that apply to it. The product meets European and national requirements on the safety of devices and products.

This product complies with the requirements of the relevant European and national directives that apply to it, RoHS.

Improper disposal of batteries/accumulators poses a risk to the environment!

Batteries/accumulators should not be disposed of with household waste. They may contain harmful heavy metals and should be treated as special waste. The chemical symbols for heavy metals are as follows: Cd = cadmium, Hg = mercury, Pb = lead.

Therefore, used batteries / accumulators should be taken to municipal hazardous waste collection points.

The built-in battery must not be disassembled for disposal. The product must be handed over in its entirety to a collection point for used electronics.

We reserve the right to make changes to the text, design and technical data of the product without notice.

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