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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

APPLICATION

The infrared massager is an advanced device designed for muscle regeneration, relaxation and improved circulation. Infrared function emits gentle heat that penetrates deep into muscle tissue, helping to relax and reduce tension. Ergonomic design and replaceable head pads allow you to adjust the massage to your individual needs.

The device is intended for home use only and is not intended for medical or commercial use and is not a substitute for medical therapy. The device is intended for use inside dry rooms. It is not intended for use in medical or therapeutic settings or for commercial use. The product is not a medical device.

For safety and CE certification reasons, the product may not be rebuilt or modified in any way. If the product is used for purposes other than those previously described, the product may be damaged. Improper use may also cause hazards such as short circuits, fire, electric shock, etc.

INSTRUCTIONS FOR USE:

Preparing the device:

- Connect the massager to a power source (AC 220-240V).

Install the appropriate cap on the head:

- Flat tip - ideal for massaging the neck, shoulders and back
- Wavy tip – contains four bulges, indispensable in the fight against fat
- Ball tip - stimulates and firms the muscle and skin
- Micro filler tip – stimulates collagen, improves skin elasticity
- Peeling tip - exfoliates the epidermis
- Tai Chi massage head - relaxes muscles, relieves tension
- Barbed tip – when used properly, it eliminates cellulite and removes excess water from the body
- Triangular roller tip – stimulates blood flow
- Three-roller tip - prepares muscles for more intensive massage techniques

Installing the overlays:

- Remove the current overlay by gently pulling it out.
- Install your chosen overlay, making sure it fits properly.

Turning on the massager:

- Turn on the device using the main switch.
- Adjust the massage intensity using the adjustment knob.

During the massage:

- Massage the selected body area using gentle, fluid movements.

After completion:

- Turn off your device.
- Unplug them from the power supply and clean the pads with a dry cloth.

TECHNICAL DATA

- Material: ABS
- Voltage: AC 220-240V, 50-60Hz
- Power: 28W
- 8 massage tips
- Adjusting the speed using the knob
- Ergonomic shape
- Infrared

SAFETY INSTRUCTIONS

- The product is intended for adults.
- Do not use the device on open wounds, swelling, bruises or highly sensitive areas.
- Do not use the device near the head or on bones.
- Keep the device out of the reach of children.
- Avoid contact of the device with water or moisture.
- Stop the massage if you experience discomfort, pain or other disturbing symptoms.
- Make sure the device is connected to a stable power source.
- Do not use the appliance if the cable or plug is damaged.
- Do not leave the device switched on unattended
- Never immerse the device in water.
- For cleaning, use a damp cloth or a mild detergent.
- Do not use a damaged device.
- The product is not intended for use by children.
- Make sure that packaging materials are not left unattended. Children may start playing with them, which is dangerous.
- Protect the product from extreme temperatures, direct sunlight, strong vibrations, high humidity, moisture, flammable gases, vapors and solvents.
- Do not expose the product to mechanical stress.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - does not function properly, - has been stored for an extended period under unfavourable conditions or - has been excessively stressed during transport.
- It is forbidden to use the product if any of its parts are damaged. In the event of damage to the cable, it is forbidden to make repairs yourself.
- Do not disassemble the device yourself.
- Protect the product from moisture.
- Turn off the product immediately and remove the cable from the product if it smells burning or smoke is visible. Have the product checked by a qualified technician before using it again.
- Before use, make sure that the existing mains voltage complies with the required operating voltage of the device.
- Connect the device only to a properly installed and grounded power outlet. The mains voltage must match the data on the nameplate of the device / mains adapter.
- Make sure that the power cable does not get wet or damp during operation. The cable must be routed so that it is not pinched or damaged.
- Keep the power cable away from hot surfaces.
- Do not perform any repairs to the device yourself.
- Any repairs may only be performed by service or by qualified specialists.
- Before cleaning the device, unplug the power supply from the wall outlet.
- Do not use the device in a damp or wet environment, such as a bathroom.
- Have a damaged plug or power cable repaired immediately by a qualified specialist or service center in order to avoid any danger.
- Do not insert any needles or sharp objects into the device.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT

Used electrical and electronic products, due to environmental protection, must not be thrown into household waste with regular municipal waste, but disposed of properly. Information on collection points and their opening hours is provided by the relevant office.

This product complies with the requirements of the relevant European and national directives that apply to it. The product meets European and national requirements on the safety of devices and products.

This product complies with the requirements of the relevant European and national directives that apply to it, RoHS.

Warning of electric shock! Danger to life!

We reserve the right to make changes to the text, design and technical data of the product without notice.

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