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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

APPLICATION AND DESCRIPTION OF THE DEVICE

The needle mesotherapy roller is ideal for performing various body treatments: on the face, belly, buttocks, thighs, etc. There are 540 stainless steel needles on the roller, which perfectly remove dead skin, while stimulating the skin to regenerate. It is also effective in the fight against cellulite. It helps in removing wrinkles, scars, stretch marks and furrows. During puncturing, the skin is stimulated to regenerate, thanks to which its condition ultimately improves and collagen production increases.

SAFETY INSTRUCTIONS

- Do not use on irritated, cut, sunburned or infected skin.
- Perform the massage gently, without excessive pressure, so as not to damage the skin.
- After treatment, disinfect the device and store it in a dry, clean place.
- Do not share the massager with other people – it is a hygienic product intended for individual use.
- If redness, burning or irritation occurs, discontinue use and consult a dermatologist.

TECHNICAL DATA

- Intensively and long-lastingly regenerates
- Improves skin elasticity, color and hydration
- Improves facial features
- Smooths scars and wrinkles
- 3 replaceable tips (180 needles - 0.50mm, 600 needles - 1.00mm, 1200 needles - 1.5mm)
- Length: 15cm
- Transparent case

HOW TO USE

- Thoroughly cleanse and dry your facial skin.
- Disinfect the roller, e.g. using isopropyl alcohol.
- Choose the appropriate attachment depending on the area you want to massage. **Smaller, narrower attachment**, designed especially for delicate areas around the eyes and mouth. **The widest roller**, designed for massaging larger parts of the face, such as cheeks, forehead, neck. **The largest roller** in the set, designed for use on the body: arms, thighs, belly.
- Gently move the massager over the skin in different directions (vertically, horizontally, diagonally).
- Don't press too hard – let the needles do the work.
- Massage of one part of the face should last about 2–3 minutes.
- When finished, apply an appropriate serum or skin care cream.
- Disinfect the device and store it in a safe, dry place.

WARNINGS

After using the device, avoid excessive sun exposure for about 1 week (it is recommended to use a cream with a high sun protection factor, e.g. 50+). Before and after each use of the device, it should be disinfected and stored in a sterile place.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center. Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the local or city office. This product complies with the requirements of the relevant European and national directives that apply to it. The product meets European and national requirements on the safety of devices and products.

We reserve the right to make changes to the text, design and technical data of the product without notice.

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