INSTRUKCJA

OBSŁUGI

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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read carefully with the instructions below for proper use of the product.

Please keep this manual for future reference and following its recommendations, because failure to follow it may pose a threat to life or health.

Introduction to the user manual

This user manual has been prepared to help you familiarize yourself with the proper and safe use of the device. It contains important instructions on operation, maintenance, and safety rules, which, if followed, will ensure the smooth operation of the device and extend its service life.

Before first use, please read the instructions carefully. The information contained herein will also help you properly store the device, especially if it will not be used for an extended period.

Following the recommendations below will allow you to enjoy trouble-free and safe operation of the product.

ATTENTION!

Using the device in a manner inconsistent with the information contained in the manual may result in the risk of bodily injury.

Application and description of the device

The neck massager is designed for at-home massage and relaxation of the neck and nape area using low-frequency electrical impulses. The device features contact electrodes that adhere to the skin and allows for a selection of different programs and intensity levels.

This device is intended for home use only and is not intended for medical or commercial use, and is not intended to replace medical therapy. The device is intended for indoor use in dry rooms. It is not intended for medical, therapeutic, or commercial use.

For safety and CE certification reasons, this product must not be rebuilt or modified in any way. Using the product for purposes other than those described may result in damage to the product. Improper use may also result in hazards such as short circuits, fire, electric shock, etc.

WARNING!

- Danger of improper use!
- Using the device in a manner other than its intended purpose and/or using the device in any other way may involve various risks.
- The device should only be used for its intended purpose.

- Please follow the procedures described in this instruction manual.
- Keep the product away from children.
- Stop the massage if you experience discomfort, pain or dizziness.

TIP

- The device must be checked for completeness of delivery and any visible damage.
- In the event of an incomplete delivery or damage due to faulty packaging or transport, please contact the service hotline.

Technical data of the device

Symbol	24128
EAN	5907451335246
Material	plastic
Power supply	battery-powered
Power	1 In
Landing	USB cable
Number of programs	12
Number of intensity levels	15
Pulse frequency	1–1000 Hz
Pulse width	20–400 µs

Startup/Use

- Please remove the packaging materials and dispose of them properly and clean the product and accessories.
- Before first use, charge the device using the USB cable. Connect the USB cable to a USB power source and connect the USB cable plug to the product's USB input port.
- Please make sure your device is properly charged before use.
- The skin in contact with the electrodes should be clean, dry and free from damage .
- Remove any jewelry from the immediate neck area.
- Gently spread the arms of the device.
- Place the massager on your neck so that both electrodes adhere evenly to the skin .
- Do not wear the device over your clothing – the electrodes should be in direct contact with the skin.
- Adjust the massager so that it does not cause pressure or discomfort.
- Turn on the device using the control button.
- Begin use at the lowest intensity level .
- Use the "+" button to gradually increase the intensity.
- Use the "-" button to decrease the intensity.
- Choose one of the 12 available work programs .
- The device allows you to set one of 15 intensity levels .
- The intensity should be increased gradually until a clear but comfortable sensation of impulses is achieved.

-
- The action should not cause sharp pain, burning or severe cramping.
 - Sit or stand in a stable position when using.
 - Do not make sudden head movements while the device is operating.
 - If the electrodes lose proper contact with the skin, turn the device off and reposition it.
 - Don't increase the intensity just to get the strongest feeling possible.
 - Once your session is over, reduce the intensity.
 - Turn off the device.
 - Only after turning it off should you remove the neck massager.
 - After use, wipe the electrodes according to the cleaning instructions.
 - Store the device in a dry place.

Intensity selection

The massager has 15 intensity levels .

- Always start at the lowest level for first use.
- Increase the intensity in single steps.
- Normal sensations may include a slight tingling or rhythmic pulses.
- The intensity should remain comfortable.
- If pain, burning, severe muscle twitching or other unpleasant sensations occur, reduce the level.
- If discomfort persists, turn off the device immediately.
- Each program may feel different, so re-evaluate the appropriate intensity level after changing modes.

Cleaning and maintenance

- Turn off the device before cleaning.
- Disconnect the USB cable.
- Do not clean the massager while charging.
- Wipe the casing with a soft, slightly damp cloth.
- Wipe the contact electrodes gently after each use.
- Do not immerse the device in water.
- Do not rinse the massager under running water.
- Do not allow water to enter the USB port.
- Do not use alcohol, gasoline, solvents or aggressive detergents.
- Do not use metal brushes or abrasive materials.
- After cleaning, allow the device to dry completely.
- Do not connect the charger while the device is wet.
- Regularly check the condition of the electrodes, housing and USB cable.
- Do not use the device with a cracked housing or damaged electrodes.
- Store the massager in a dry, clean place.

Safety tips and rules

- Please read the entire manual before first use.
- Use the device only for its intended purpose.
- The massager is intended for personal use.
- Do not use the device as a substitute for treatment or professional medical equipment.
- Do not use the device on the front of the neck or throat area.
- The electrodes should only be placed in the position intended by the design on the back/side of the neck.
- Do not use the device on the head, chest or directly in the heart area.
- Do not use the massager on skin with wounds, abrasions, burns, inflammation or other visible damage.
- Do not use it on areas with significantly reduced sensation.
- Do not use the device simultaneously with another device that uses electrical impulses.
- People with a pacemaker, implanted defibrillator or other active electronic implant should not use the device without consulting their physician .
- People with heart disease or heart rhythm disorders should consult a doctor before use.
- Pregnant women should consult their doctor before using the device.
- If you have any doubts about your health, consult a specialist before use.
- Do not use the device while driving.
- Do not use it while riding a bicycle or operating machinery.
- Do not use the device while sleeping.
- Don't fall asleep with the massager on.
- Do not use the device in the bathtub, shower or other places with high humidity.
- Do not use the device with wet hands.
- Do not immerse the device in water.
- Do not expose the device to rain.
- Do not use the massager while charging.
- Do not connect a damaged USB cable.
- Do not insert any objects into the charging socket.
- Protect the device from high temperatures and open flames.
- Do not leave the massager on a radiator or other heat source.
- Do not expose the device to direct sunlight for long periods of time.
- Do not drop the device.
- Do not hit or crush its casing.
- Do not bend the arms of the device with excessive force.
- Do not attempt to disassemble the device yourself.
- Do not modify the electronics or electrodes.
- Before putting on the device, check the condition of your skin.
- The electrodes must adhere evenly to the skin.
- Do not increase the intensity suddenly.
- Always start with minimum intensity when using for the first time.

- The impulses should not cause sharp pain.
- If you experience burning, severe stinging, dizziness, feeling unwell or other disturbing symptoms, turn off the device immediately.
- In case of persistent discomfort, discontinue use.
- Do not use the device in a way that causes severe, painful muscle contractions.
- Do not wear the device over thick clothing.
- Do not place metal objects between the electrodes and your skin.
- Before use, remove any metal jewelry that may be in the area of contact.
- Do not use the device immediately after applying flammable or alcohol-based preparations to the skin.
- Keep the device out of the reach of small children.
- Children should not use the device as a toy.
- Do not allow children to regulate the intensity of impulses on their own.
- Do not use the device if it is damaged or operates abnormally.
- After use, turn off the device and store it in a dry, safe place.

Safety instructions for batteries/rechargeable batteries

- Keep batteries/rechargeable batteries out of the reach of children. If swallowed, seek medical attention immediately!
- Disposable batteries must not be recharged. Batteries/rechargeable batteries should not be short-circuited and/or opened. This may result in overheating, fire, or explosion.
- Never throw batteries/rechargeable batteries into fire or water.
- Batteries/rechargeable batteries should never be subjected to mechanical stress.
- Risk of acid leakage from batteries/accumulators.
- Avoid extreme conditions and temperatures that may affect batteries/accumulators, e.g. radiators/direct sunlight.
- If batteries/rechargeable batteries leak, avoid contact with skin, eyes, and mucous membranes! Immediately flush the affected areas with clean water and seek medical attention!
- Leaking or damaged batteries/rechargeable batteries can cause chemical burns when in contact with skin. Therefore, wear appropriate protective gloves.
- The product has a built-in battery that cannot be replaced by the user.

Warning and information symbols

	<u>TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING</u> The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center. Used packaging material should be delivered to a designated waste disposal site designated by local authorities. For information on recycling options for used products, please contact your local municipality or city office.
---	---

	<u>DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT</u> For environmental reasons, used electrical and electronic products should not be disposed of with regular household waste, but should be disposed of properly. Information on collection points and their opening hours can be obtained from the relevant authorities.
	This product complies with the requirements of the applicable European and national directives. The product meets European and national requirements for the safety of devices and products.
	This product complies with the requirements of the relevant European and national RoHS directives.
	<u>Improper disposal of batteries/accumulators poses a threat to the environment!</u> Batteries/accumulators should not be disposed of with household waste. They may contain harmful heavy metals and should be treated as special waste. The chemical symbols for heavy metals are: Cd = cadmium, Hg = mercury, Pb = lead. Therefore, used batteries/accumulators should be taken to municipal hazardous waste collection points. The built-in battery must not be disassembled for disposal. The product should be taken to a collection point for used electronics.

We reserve the right to make changes to the text, design and product data without notice.

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