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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

APPLICATION

Thanks to the technology of deep percussion massage, the gun effectively reduces muscle tension, improves circulation and supports regeneration. Various tips allow you to adjust the massage to specific parts of the body, which ensures comprehensive muscle care. Professional device for massaging muscles, back, legs and the whole body! The massage gun is an ideal solution for physically active people, athletes, as well as all those who feel muscle tension and fatigue. Thanks to the multifunctional design and 8 interchangeable tips, it allows you to adjust the massage to different parts of the body. The massage gun is mainly used to relieve muscle stiffness and pain, to improve blood circulation. It improves the condition of soft tissues, accelerates and regenerates inflammation. It is perfect for people who engage in recreational physical activity as well as for professional athletes. It works best after a hard workout or after long physical work.

INSTRUCTIONS FOR USE:

Charge your device using the included USB cable.

Choose the appropriate tip and place it in the gun socket. (ROUND: universal - designed for massaging all parts of the body, e.g. arms, legs, back, buttocks. FLAT: designed for working on small and large muscle groups, relaxing. PROJECTILE: designed for trigger points (thickenings on the tissue that cause pain), for deep tissue, joints. FORK: designed for working on narrower and long muscles (e.g. erector spinae, trapezius, gastrocnemius, Achilles).

Turn on the device by pressing the power button.

Choose the appropriate intensity level, adapting the massage to your needs.

Place the massager on the selected body part (arms, back, thighs, calves) and slowly move it along the muscles. Do not press the device too hard - its weight and movement are enough for an effective massage.

One massage session for a given body part should last a maximum of 10–15 minutes. The total time of using the device during the day should not exceed 30 minutes.

Once you have finished massaging, turn the device off by holding the power button.

Clean the massage tip with a dry, soft cloth and store the device in a dry place.

TIP

- The device must be checked for completeness of delivery and any visible damage.
- In the event of an incomplete delivery or if you notice damage due to faulty packaging or transport, please contact the service hotline.

TECHNICAL DATA

- Touch control panel
- 10 speed levels
- One work cycle lasts 15 minutes
- Battery: lithium-ion
- Charging time about 120 minutes
- Rotation speed: 2600 – 3200r/min
- Power: 24W
- 4 sets of different tips
- Power supply: built-in battery
- Charging: USB cable – Type-C, 5V/1A
- Dimensions without tips: 18cm x 14.5cm x 10cm

SAFETY INSTRUCTIONS

- This product may be used by children aged 8 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- Do not use the device on open wounds, swelling, bruises or sensitive areas.
- Do not use the device on the head, neck or bones.
- Keep the device out of the reach of children.
- Avoid contact of the device with water or moisture.
- Stop the massage if you experience discomfort, pain or other disturbing symptoms.
- Only charge your device using the original USB cable.
- Do not leave the device switched on unattended
- Never immerse the device in water.
- For cleaning, use a damp cloth or a mild detergent.
- Do not use a damaged device.
- The product is not intended for use by children.

- Make sure that packaging materials are not left unattended. Children may start playing with them, which is dangerous.
- Protect the product from extreme temperatures, direct sunlight, strong vibrations, high humidity, moisture, flammable gases, vapors and solvents.
- Do not expose the product to mechanical stress.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - does not function properly, - has been stored for an extended period under unfavourable conditions or - has been excessively stressed during transport.
- It is forbidden to use the product if any of its parts are damaged. In the event of damage to the cable, it is forbidden to make repairs yourself.
- Do not disassemble the device yourself.
- Protect the product from moisture.
- Make sure not to damage the USB cable with sharp edges or hot objects. Fully unwind the USB cable before use.
- The included USB cable is only suitable for this product.
- Immediately turn off the product and remove the cable from the product if it smells burning or smoke is visible. Have the product checked by a qualified technician before using it again.
- If a USB power supply is used, the socket must always be easily accessible so that the USB power supply can be easily removed from the socket in an emergency. Please also observe the manufacturer's operating instructions.
- Before each charging process, check the USB cable and device for any damage. Never use or charge the device if you find any damage.
- Before use, make sure that the existing mains voltage matches the required operating voltage of the device.

BATTERIES/ACCUMULATORS SAFETY INSTRUCTIONS

- Batteries/rechargeable batteries should be kept out of reach of children. If swallowed, consult a doctor immediately!
- Disposable batteries must not be recharged. Batteries/accumulators must not be short-circuited and/or opened. This may cause overheating, fire or explosion.
- Never throw batteries/accumulators into fire or water.
- Batteries/accumulators should never be subjected to mechanical stress.
- Risk of acid leakage from batteries/accumulators.
- Avoid extreme conditions and temperatures that may affect batteries / accumulators, e.g. radiators / direct sunlight.
- If batteries/accumulators leak, avoid contact of skin, eyes and mucous membranes with chemicals! Immediately rinse affected areas with clean water and seek medical attention!
- Leaking or damaged batteries/rechargeable batteries can cause chemical burns when in contact with skin. Therefore, suitable protective gloves should be worn in such cases.
- The product has a built-in battery that cannot be replaced by the user.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT

Used electrical and electronic products, due to environmental protection, must not be thrown into household waste with regular municipal waste, but disposed of properly. Information on collection points and their opening hours is provided by the appropriate office.

This product complies with the requirements of the relevant European and national directives that apply to it. The product meets European and national requirements on the safety of devices and products.

This product complies with the requirements of the relevant European and national directives that apply to it, RoHS.

Improper disposal of batteries/accumulators poses a risk to the environment!

Batteries/accumulators should not be disposed of with household waste. They may contain harmful heavy metals and should be treated as special waste. The chemical symbols for heavy metals are as follows: Cd = cadmium, Hg = mercury, Pb = lead.

Therefore, used batteries / accumulators should be taken to municipal hazardous waste collection points. The built-in battery must not be disassembled for disposal. The product must be handed over in its entirety to a collection point for used electronics.

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