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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow them may pose a risk to life or health.

APPLICATION AND DESCRIPTION

The Roller Foot Massager is a device designed for a relaxing and stimulating massage of the soles of the feet. It features two sets of massage rollers and metal balls placed at the ends of the massage modules. The movement of the rollers stimulates the receptors in the feet, improves microcirculation, relieves tension, and supports overall body regeneration. Its compact size makes it easy to use at home, in the office, or while traveling.

TIP

- The device must be checked for completeness of delivery and any visible damage.
- In the event of an incomplete delivery or damage due to faulty packaging or transport, please contact the service hotline.

TECHNICAL DATA

- Product material: plastic, metal
- Massage rollers: 16 pieces
- Massage type: acupressure, rotational massage

SAFETY INSTRUCTIONS

- Before installing and using the product, please read this instruction manual carefully.
- This product may be used by children aged 8 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- Children must not play with the product.
- Cleaning and maintenance must not be performed by children without supervision.
- Make sure packaging materials are not left unattended. Children may play with them, which is dangerous.
- Keep the foil packaging away from children (risk of suffocation)
- Avoid sources of ignition: avoid contact with fire or heat-generating devices (may cause fire)
- Protect the product from extreme temperatures, direct sunlight, high humidity, flammable gases, and solvents.
- For cleaning, use a damp cloth or a mild detergent.
- Do not use the device on damaged, irritated or injured skin on the feet.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - is not functioning properly, - has been stored for an extended period under unfavorable conditions, or - has been subjected to excessive stress during transport.
- If you experience severe pain, numbness or discomfort, stop the massage.
- People with circulatory problems, diabetes or foot conditions should consult a physician before use.
- Do not stand with your full body weight on the device – the massager is not a platform for standing.
- The product should always be used as intended.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a designated waste disposal site designated by local authorities. For information on recycling options for used products, please contact your local municipality or city office.

The product meets European and national requirements on the safety of devices and products.
We reserve the right to make changes to the text, design and product data without notice.

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