- (PL) INSTRUKCJA OBSŁUGI - WSKAZÓWKI DOTYCZĄCE OBSŁUGI I BEZPIECZEŃSTWA**
OPASKA NA KARK NY-93 SYMBOL: 27289 EAN/GTIN:5907451364598
- (EN) USER MANUAL - OPERATING AND SAFETY INSTRUCTIONS**
NECKBAND SYMBOL: 27289 EAN/GTIN:5907451364598
- (DE) BEDIENUNGSANLEITUNG - BEDIENUNGS- UND SICHERHEITSHINWEISE**
NACKENBAND SYMBOL: 27289 EAN/GTIN:5907451364598
- (CZ) NÁVOD K OBSLUZE - PROVOZNÍ A BEZPEČNOSTNÍ RADY**
NÁKRČNÍ PÁSEK SYMBOL: 27289 EAN/GTIN:5907451364598
- (FR) MANUEL D'UTILISATION - CONSEILS D'UTILISATION ET DE SÉCURITÉ**
BANDEAU DE NUQUE SYMBOLE: 27289 EAN/GTIN:5907451364598
- (IT) MANUALE OPERATIVO - CONSIGLI OPERATIVI E DI SICUREZZA**
FASCIA PER IL COLLO SIMBOLO: 27289 EAN/GTIN:5907451364598
- (ES) MANUAL DE OPERACIÓN - CONSEJOS DE OPERACIÓN Y SEGURIDAD**
BANDA PARA EL CUELLO SÍMBOLO: 27289 EAN/GTIN:5907451364598
- (NL) BEDIENINGSHANDLEIDING - BEDIENINGS- EN VEILIGHEIDTIPS**
NEKBAND SYMBOOL: 27289 EAN/GTIN:5907451364598
- (SE) BRUKSANVISNING - DRIFT- OCH SÄKERHETSTIPS**
NACKBAND SYMBOL: 27289 EAN/GTIN:5907451364598
- (GR) ΕΓΧΕΙΡΙΔΙΟ ΛΕΙΤΟΥΡΓΙΑΣ - ΣΥΜΒΟΥΛΕΣ ΛΕΙΤΟΥΡΓΙΑΣ ΚΑΙ ΑΣΦΑΛΕΙΑΣ**
ΚΟΡΔΕΛΑ ΓΙΑ ΤΟΝ ΑΥΧΕΝΑ ΣΥΜΒΟΛΟ: 27289 EAN/GTIN:5907451364598
- (RO) MANUAL DE OPERARE - SFATURI DE OPERAREA SI SIGURANTA**
BANDĂ PENTRU GÂT SIMBOL: 27289 EAN/GTIN:5907451364598
- (PT) MANUAL DE OPERAÇÃO - DICAS DE OPERAÇÃO E SEGURANÇA**
FAIXA DE PESCOÇO SÍMBOLO: 27289 EAN/GTIN:5907451364598
- (BG) РЪКОВОДСТВО ЗА ЕКСПЛОАТАЦИЯ - СЪВЕТИ ЗА ЕКСПЛОАТАЦИЯ И БЕЗОПАСНОСТ**
ЛЕХТА ЗА ВРАТА СИМБОЛ: 27289 EAN/GTIN:5907451364598
- (HU) HASZNÁLATI ÚTMUTATÓ – HASZNÁLATI ÉS BIZTONSÁGI TIPP EK**
NYAKPÁNT SZIMBÓLUM: 27289 EAN/GTIN:5907451364598
- (DK) BETJENINGSVEJLEDNING - BETJENINGS- OG SIKKERHEDSTIPS**
NAKKEBÅND SYMBOL: 27289 EAN/GTIN:5907451364598
- (SK) NÁVOD NA OBSLUHU - RADY NA PREVÁDZKU A BEZPEČNOSŤ**
PÁS NA KRK SYMBOL: 27289 EAN/GTIN:5907451364598
- (FI) KÄYTTÖOHJE - KÄYTTÖ- JA TURVALLISUUSVINKKEJÄ**
NISKAPANTA SYMBOLI: 27289 EAN/GTIN:5907451364598
- (LT) NAUDOJIMO VADOVAS – NAUDOJIMO IR SAUGOS PATARIMAI**
JUOSTA ANT KAKLO SIMBOLIS: 27289 EAN/GTIN:5907451364598
- (LV) LIETOŠANAS ROKASGRĀMATA – LIETOŠANAS UN DROŠĪBAS PADOMI**
JOSTA KAKLAM SIMBOLS: 27289 EAN/GTIN:5907451364598
- (EE) KASUTUSJUHEND – KASUTUS- JA OHUTUSJUHEND**
KAELAPAEEL SÜMBOL: 27289 EAN/GTIN:5907451364598
- (SI) PRIROČNIK ZA UPORABO - NASVETI ZA UPORABO IN VARNOST**
TRAK ZA VRAT SIMBOL: 27289 EAN/GTIN:5907451364598
- (IE) LÁMHLEABHAR OIBRIÚCHÁIN - TIPS OIBRIÚCHÁIN AGUS SÁBHÁILTEACHTA**
BANDA MUINÉIL SIOMBOOL: 27289 EAN/GTIN:5907451364598
- (MT) MANWAL OPERAT - OPERAZZJONI U GHAJR TAS-SIGURTÀ**
FAXXA GHALL-GHATX SIMBOLU: 27289 EAN/GTIN:5907451364598
- (HR) PRIRUČNIK ZA UPORABU - SAVJETI ZA UPOTREBU I SIGURNOST**
TRAKA ZA VRAT SIMBOL: 27289 EAN/GTIN:5907451364598
- (RU) РУКОВОДСТВО ПО ЭКСПЛУАТАЦИИ – СОВЕТЫ ПО ЭКСПЛУАТАЦИИ И БЕЗОПАСНОСТИ**
ПОВЯЗКА НА ШЕЮ СИМБОЛ: 27289 EAN/GTIN:5907451364598

Importer / Importeur / Importateur / Importatore / Importador / Importör / Εισαγωγέας / Importator / Dovožce / Importador / Вноситель / Importör / Importør / Dovožca / Maahantuojaja / Importuotojas / Importētājs / Maaletooja / Uvoznik / Allmhaireoir / Importatur / Uvoznik / Импортёр:

VERK GROUP SIKORSKI SP.K.

WYGODY 16

05-090 PODOLSZYN NOWY

POLAND

info@verkgroup.pl

+48 888055550 / +48 22 670 90 10

internetowa-hurtownia.pl

verk.store

verk.pl

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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

APPLICATION AND DESCRIPTION OF THE DEVICE

The neck and neck massage band is a compact and functional device that helps relieve pain, muscle tension and stress. Equipped with a warming function, it effectively supports muscle relaxation and improves blood circulation, bringing relief after a long day at work or intense exercise.

Thanks to its USB power supply, the headband is comfortable and mobile - it can be used at home, in the office or on the go. The lightweight design and ergonomic shape ensure comfort when worn for a longer period of time. The massager offers different vibration and heating strengths, which can be adjusted to individual needs, ensuring precise action in the critical areas of the neck and nape.

This practical device combines relaxation with regeneration, ideal for those struggling with muscle pain or looking for a daily moment of respite. The massage band is easy to use, compact and stylish, making it a great choice for your own use or as a gift for a loved one.

TIP

The device must be checked for completeness of delivery and any visible damage.

In the event of an incomplete delivery or if you notice damage due to faulty packaging or transport, please contact the service hotline.

TECHNICAL DATA

- Product material: ABS material
- Bluetooth support
- Controlled also by remote control
- Product size: ~ 12 x 16 x 16cm
- 3 brightness levels
- Power supply: USB DC 5V
- Height: 14cm
- Strap length: 26cm
- Total length: 65cm

SAFETY INSTRUCTIONS

- This product may be used by children aged 8 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- Children are not allowed to play with the product.
- Cleaning and maintenance may not be performed by children without supervision.
- Never immerse the device in water.
- For cleaning, use a damp cloth or a mild detergent.
- Do not use a damaged device.
- Make sure that packaging materials are not left unattended. Children may start playing with them, which is dangerous.
- Protect the product from extreme temperatures, direct sunlight, strong vibrations, high humidity, moisture, flammable gases, vapors and solvents.
- Do not expose the product to mechanical stress.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - does not function properly, - has been stored for an extended period under unfavourable conditions or - has been excessively stressed during transport.
- It is forbidden to use the product if any of its parts are damaged. In the event of damage to the cable, it is forbidden to make repairs yourself.
- Do not disassemble the device yourself.
- Protect the product from moisture.
- The light elements cannot be replaced. If the light elements have stopped working due to wear, the entire product must be replaced.
- Make sure not to damage the USB cable with sharp edges or hot objects. Fully unwind the USB cable before use.
- The included USB cable is only suitable for this product.
- Immediately turn off the product and remove the cable from the product if it smells burning or smoke is visible. Have the product checked by a qualified technician before using it again.

- If a USB power supply is used, the socket must always be easily accessible so that the USB power supply can be easily removed from the socket in an emergency. Please also observe the manufacturer's operating instructions.
- Before use, make sure that the existing mains voltage matches the required operating voltage of the device.
- Disconnect the USB cable before assembly, disassembly and cleaning.
- Always use the product as intended.

BATTERIES/ACCUMULATORS SAFETY INSTRUCTIONS

- Batteries/rechargeable batteries should be kept out of reach of children. If swallowed, consult a doctor immediately!
- Disposable batteries must not be recharged. Batteries/accumulators must not be short-circuited and/or opened. This may cause overheating, fire or explosion.
- Never throw batteries/accumulators into fire or water.
- Batteries/accumulators should never be subjected to mechanical stress.
- Risk of acid leakage from batteries/accumulators.
- Avoid extreme conditions and temperatures that may affect batteries / accumulators, e.g. radiators / direct sunlight.
- If batteries/accumulators leak, avoid contact of skin, eyes and mucous membranes with chemicals! Immediately rinse affected areas with clean water and seek medical attention!
- Leaking or damaged batteries/rechargeable batteries can cause chemical burns when in contact with skin. Therefore, suitable protective gloves should be worn in such cases.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT

Used electrical and electronic products, due to environmental protection, must not be thrown into household waste with regular municipal waste, but disposed of properly. Information on collection points and their opening hours is provided by the relevant office.

This product complies with the requirements of the relevant European and national directives that apply to it. The product meets European and national requirements on the safety of devices and products.

This product complies with the requirements of the relevant European and national directives that apply to it, RoHS.

Do not look into a light source when it is on. It may be harmful to your eyes.

Improper disposal of batteries/accumulators poses a risk to the environment!

Batteries/accumulators should not be disposed of with household waste. They may contain harmful heavy metals and should be treated as special waste. The chemical symbols for heavy metals are as follows: Cd = cadmium, Hg = mercury, Pb = lead.

Therefore, used batteries / accumulators should be taken to municipal hazardous waste collection points.

The built-in battery must not be disassembled for disposal. The product must be handed over in its entirety to a collection point for used electronics.

We reserve the right to make changes to the text, design and technical data of the product without notice.

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